

# **FID-A: The FID Appliance**

Version 1.0

User Manual

Jamie Near

10 December 2014

## TABLE OF CONTENTS

TABLE OF CONTENTS	2
1. INTRODUCTION	7
1.1. Simulation Toolbox	7
1.2. RF Pulse Toolbox	8
1.3. Processing Toolbox	9
1.4. Data structure formatting	10
2. SIMULATION TOOLS	13
2.1. sim_Hamiltonian.m	13
2.2. sim_evolve.m	13
2.3. sim_excite.m	14
2.4. sim_excite_arbPh.m	14
2.5. sim_lcmrawbasis.m	14
2.6. sim_megapress.m	15
2.7. sim_megapress_shaped.m	16
2.8. sim_megapress_shapedEdit.m	17
2.9. sim_megapress_shapedRefoc.m	17
2.10. sim_onepulse.m	18
2.11. sim_onepulse_arbPh.m	19
2.12. sim_onepulse_shaped.m	19
2.13. sim_press.m	20

2.14. sim_press_shaped.m	20
2.15. sim_readout.m	21
2.16. sim_rotate.m	21
2.17. sim_rotate_arbPh.m	22
2.18. sim_runAllMetabs.m	22
2.19. sim_runMegaPressShaped.m	22
2.20. sim_runMegaPressShapedEdit.m	23
2.21. sim_runMegaPressShapedRefoc.m	24
2.22. sim_shapedRF.m	25
2.23. sim_spinecho.m	25
2.24. sim_spinecho_shaped.m	25
2.25. sim_spoil.m	26
2.26. sim_steam.m	26
2.27. sim_steam_gradSim.m	26
2.28. writelcmraw.m	27
<b>3. RF PULSE TOOLS</b>	<b>27</b>
3.1. rf_blochSim.m	27
3.2. rf_dualBand.m	27
3.3. rf_freqshift.m	28
3.4. rf_gauss.m	28
3.5. rf_hs.m	28
3.6. rf_loadwaveform.m	29
3.7. rf_readRF.m	29
3.8. rf_readpta.m	29
3.9. rf_resample.m	30
3.10. rf_writepta.m	30
3.11. rf_writerf.m	30
<b>4. PROCESSING TOOLS</b>	<b>30</b>
4.1. addphase.m	30

4.2. addphase1.m	31
4.3. getLW.m	31
4.4. getLWandSNR.m	31
4.5. getSNR.m	32
4.6. op_ISIScombine.m	32
4.7. op_addNoise.m	32
4.8. op_addScans.m	33
4.9. op_addphase.m	33
4.10. op_addphaseSubspec.m	33
4.11. op_addrcvrs.m	33
4.12. op_alignAverages.m	34
4.13. op_alignAverages_fd.m	34
4.14. op_alignScans.m	34
4.15. op_alignScans_fd.m	35
4.16. op_alignrcvrs.m	35
4.17. op_ampScale.m	35
4.18. op_averaging.m	36
4.19. op_avgNormalize.m	36
4.20. op_combinesubspecs.m	36
4.21. op_complexConj.m	37
4.22. op_concatAverages.m	37
4.23. op_creFit.m	37
4.24. op_dccorr.m	37
4.25. op_downsamp.m	38
4.26. op_ecc.m	38
4.27. op_fddccorr.m	38
4.28. op_filter.m	38
4.29. op_freqAlignAverages.m	39
4.30. op_freqAlignAverages_fd.m	39

4.31. op_freqrange.m	39
4.32. op_freqshift.m	40
4.33. op_freqshiftSubspec.m	40
4.34. op_gaussianPeak.m	40
4.35. op_getcoilcombos.m	40
4.36. op_getcoilcombos_specReg.m	41
4.37. op_leftshift.m	41
4.38. op_loadjmrui.m	41
4.39. op_loadlcmdetail.m	42
4.40. op_loadspec_IMA.m	42
4.41. op_loadspec_twix.m	42
4.42. op_loadspec_varian.m	42
4.43. op_lorentz.m	43
4.44. op_lorentz_linbas.m	43
4.45. op_lorentzianPeak.m	44
4.46. op_makeFreqDrift.m	44
4.47. op_makePhaseDrift.m	44
4.48. op_megapressproc.m	44
4.49. op_phaseAlignAverages.m	45
4.50. op_phaseAlignAverages_fd.m	45
4.51. op_pressproc.m	46
4.52. op_readlcmraw.m	46
4.53. op_readlcmraw_basis.m	46
4.54. op_rmNworstaverages.m	47
4.55. op_rmbadaverages.m	47
4.56. op_rmworstaverage.m	47
4.57. op_specialproc.m	47
4.58. op_specialproc_fmrs.m	48
4.59. op_specialproc_fmrs_slidingWindow.m	49

4.60. op_subtractScans.m	49
4.61. op_takeaverages.m	49
4.62. op_takesubspec.m	50
4.63. op_unfilter.m	50
4.64. op_writejmrui.m	50
4.65. op_writelcm.m	50
4.66. op_zeropad.m	51
4.67. op_zerotrim.m	51
4.68. readlcmcoord.m	51
4.69. readlcmcoord_getBackground.m	51
4.70. readlcmtab.m	52

# 1. Introduction

The FID Appliance (FID-A) is an open-source software package for simulation of MRS experiments, design and analysis of radiofrequency (RF) pulses, and processing of MRS data. The software is freely available for download ([www.github.com/CIC-methods/FID-A](http://www.github.com/CIC-methods/FID-A)).

The FID-A software package consists of three separate toolboxes; the Simulation Toolbox, the RF-pulse toolbox, and the Processing toolbox. Each is described briefly below:

## 1.1. Simulation Toolbox

The FID-A NMR simulation toolbox is based on an implementation of the density matrix formalism, where the evolution of the spin system in a given NMR experiment is described by successive evolutions of the density matrix by time-independent Hamiltonian operators.

The simulation toolbox contains built-in functions for simulating the basic components of an MRS pulse sequence, including excitation (`sim_excite.m`), delays (`sim_evolve.m`), rotation using ideal (`sim_rotate.m`) and shaped (`sim_shapedRF.m`) radiofrequency pulses, and signal readout (`sim_readout.m`). Using these basic pulse sequence elements, common in-vivo MRS pulse sequences are also implemented, including the FID (`sim_onepulse.m`), PRESS (`sim_press.m`), STEAM (`sim_steam.m`) and spin-echo (`sim_spinecho.m`).

Simulation is performed on a given spin-system by specifying the pulse timings and the chemical shifts and coupling constants of the spin-system of interest. A full set of common metabolite spin-system definitions is provided based on values previously

published by Govindaraju et al (NMR Biomed 2000, 13:129-153). Spin system definitions can be found in .../FID-A/simulationTools/metabolites/. Each metabolite's spin system definition is contained within a structure with two fields: "shifts" and "J". Shifts is a vector of length N, which contains the chemical shift value (in ppm) for each of the N spins in the spin system. J is an N x N matrix specifying the coupling constant (in Hz) between each spin and every other spin in the system. For spin systems with more than 7 spins, the spin system is broken down in to multiple independent (uncoupled) components with fewer spins. This is because the computation time increases more than linearly with increasing number of spins, and so it is much more computationally efficient to simulate, say two groups of spins with 7 protons each rather than simulating one group of 14 protons. Finally, some spin systems have been reduced in size due to redundancy. For example, scyllo inositol consists of nine magnetically equivalent protons, and is therefore most easily simulated as a single proton. However, the user must then remember to scale the intensity of Scyllo Inositol Simulations by a factor of 9 using `op_ampscale(sim_scyllo,9);`.

Excitation and refocusing RF pulses can be modeled as ideal (instantaneous) rotations, or fully shaped RF waveforms, depending on the user requirements. In the case of shaped RF pulses, phase cycling is performed to remove unwanted coherences. Function names in the simulation toolbox begin with the prefix "sim\_".

## **1.2. RF Pulse Toolbox**

The FID-A RF toolbox enables the creation of basic RF pulse waveforms, and Bloch simulation to determine the resulting excitation/refocusing/inversion profiles, as well as frequency



shifting and resampling of rf waveforms. RF pulses are stored in MATLAB structure format, with fields corresponding to the RF waveform, the type of rf pulse (excitation, inversion, refocusing), the time-bandwidth product of the rf pulse, and the time-B1 product. Function names in the RF-pulse toolbox begin with the prefix "rf\_".

### **1.3. Processing Toolbox**

The FID-A processing toolbox contains "load" functions to accept MRS data in MRI vendor data formats (Siemens, Agilent, others in development), and store them in MATLAB.

The data can then be operated on using any of the over 50 different processing operations, including (but not limited to) filtering (op\_filter.m), zeropadding (op\_zeropad.m), time domain truncation (op\_leftshift.m and op\_zerotrim.m), frequency domain truncation (op\_freqrange.m), eddy current correction (op\_ecc.m), removal motion corrupted averages (op\_rmbadaverages.m), retrospective frequency and phase drift correction (op\_alignAverages.m and op\_alignAverages\_fd.m, Near et al, Magn Reson Med 2014, DOI: 10.1002/mrm.25094), combination of multi-element RF coil data (op\_addrcvrs.m), and zero- and first-order phase corrections (op\_addphase.m).

These functions can be nested within one another. For example, to load a spectrum, combine the receivers, combine the averages and then filter the result, can be done using the following single line of code:

```
out=
op_filter(op_averaging(op_addrcvrs(op_loadspec_twix('filename.dat'),1,'w'),5);
```

where the argument "5" represents a 5 Hz exponential filter, the argument "'w'" specifies that a weighted coil recombination should be performed, and the argument "1" specifies that the phases and amplitudes of the rf coil channels should be determined using the first point in the time domain.

Finally, the FID-A processing toolbox contains "write" functions to output the data into the formats accepted by leading MRS data processing packages. Function names in the processing toolbox begin with the prefix "op\_", which stands for operator.

While most of the tools in the processing toolbox are single processing operations, there are a few examples of full MRS data processing pipelines. These include: `op_specialproc.m`, `op_pressproc.m`, and `op_megapressproc.m`. These "pipeline" operations begin with raw MRS data in Siemens .dat format and then process the data in a logical step-by-step fashion to generate fully processed data that is ready to be analyzed using one of the leading MRS analysis software packages (LCModel, jMRUI or Tarquin). Processing steps performed in these pipelines include combination of receive channels (`op_addrcvrs.m`), removal of motion corrupted averages (`op_rmbadaverages.m`), spectral registration of averages (`op_alignAverages.m`), and combination of averages (`op_averaging.m`). For more information on these pipeline operations, see section 4.

#### **1.4. Data structure formatting**

The FID-A software package is implemented in MATLAB (Natick MA, USA). Within FID-A, simulated or experimental MRS datasets, and RF pulses are stored in uniquely formatted data structures that encapsulate all of the data and relevant header information. Each dataset structure stores both time-domain and frequency-domain data arrays. When multiple averages, coil channels, or

subspectra of data are present, they are stored in separate dimensions of the arrays, and indexed within the header. Each RF pulse structure stores the rf waveform as well as information about the rf pulse type (excitation, refocusing or inversion), the time-bandwidth product, and the time-B1 product. By encapsulating data and header information in this way, processing operations have access to all relevant information and thus require as few input arguments as possible. The fields of the data structure are listed and briefly described below:

fids	- time domain MRS data.
specs	- frequency domain MRS data.
t	- vector of time values for plotting in the time domain [s]
ppm	- vector of frequency values for plotting in the frequency domain [ppm]
sz	- size of the fids and specs arrays
date	- date that the data was acquired or simulated
averages	- number of averages in the dataset (possibly altered by processing)
rawAverages	- number of averages in the original dataset (not altered by processing).
subspecs	- number of subspectra (ISIS, edit on/off, etc) in the dataset (possibly altered by processing).
rawSubspecs	- number of subspectra (ISIS, edit on/off, etc) in the original dataset (not altered by processing).
Bo	- magnetic field strength [Tesla]
txfrq	- Centre frequency [MHz];
linewidth	- linewidth of data (only used for simulated data) [Hz]
n	- number of spectral points
dweltime	- dwell time of the data in the time domain [s] (dweltime = 1/spectralwidth)
sim	- type of simulation (ideal vs. shaped pulses), only used for simulated data.
te	- echo time of acquisition [ms], only used for simulated data
seq	- type of sequence used (only used for simulated data).
dims	- structure specifying which data dimensions are stored along which dimensions of the fids/specs arrays. Fields include:
t	- time/frequency dimension (usually this is 1, the first dimension of the fids/specs array).



the rf pulse [Hz].  
tbw        - The product of the duration of the rf pulse [s] and the  
            bandwidth of the rf pulse [Hz]

Below is a brief summary of all of the functions in the FID-A  
Software package. The same information can be obtained by typing  
'help functionName' at the MATLAB command line:

## **2. Simulation Tools**

### **2.1. sim\_Hamiltonian.m**

**USAGE:**

[H,d] = sim\_Hamiltonian(sys,Bfield);

**DESCRIPTION:**

Creates the nxn Hamiltonian matrix for a spin system, which can then be used  
in other functions to simulate NMR experiments.

**INPUTS:**

sys                = spin system definition structure.  
Bfield            = magnetic field strength (Tesla).

### **2.2. sim\_evolve.m**

**USAGE:**

d\_out = sim\_evolve(d\_in,H,t)

**DESCRIPTION:**

This function simulates free evolution of the spin system under the effects of  
chemical shift and scalar coupling.

**INPUTS:**

d\_in            = input density matrix structure.  
H               = Hamiltonian operator structure.  
t               = duration of evolution (s)

### 2.3. **sim\_excite.m**

**USAGE:**

```
d_out = sim_excite(H,axis,angle)
```

**DESCRIPTION:**

This function simulates the effect of an ideal (instantaneous) excitation pulse on the density matrix. Used in simulation tools.

**INPUTS:**

```
d_in      = input density matrix structure.  
H         = Hamiltonian operator structure.  
axis      = Axis of rotation ('x' or 'y');  
angle     = Flip angle of excitation (degrees). Optional. Default=90. If  
            angle is a scalar, then the same flip angle is applied to all  
            spins in the spin system. If angle is a vector, then the elements  
            of the vector specify the flip angles to apply to each spin in the  
            system. In this case, the length of the vector must be the same  
            as the number of spins in the spin system.
```

### 2.4. **sim\_excite\_arbPh.m**

**USAGE:**

```
d_out = sim_excite_arbPh(H,phase,angle)
```

**DESCRIPTION:**

This function simulates the effect of an ideal (instantaneous) excitation pulse on the density matrix. Used in simulation tools. The phase of the excitation pulse can be arbitrarily chosen. To achieve an arbitrary phase, this code executes a rotation about z by an angle of -phi (the rf pulse phase), then executes a rotation about x by an angle of "angle" (the flip angle of the excitation pulse), and then executes a rotation back about z by an angle of phi.

**INPUTS:**

```
d_in      = input density matrix structure.  
H         = Hamiltonian operator structure.  
phase     = Phase of rotation in degrees (ie. 0='x', 90='y', etc);  
angle     = Flip angle of excitation (degrees). Optional. Default=90. If  
            angle is a scalar, then the same flip angle is applied to all  
            spins in the spin system. If angle is a vector, then the elements  
            of the vector specify the flip angles to apply to each spin in the  
            system. In this case, the length of the vector must be the same  
            as the number of spins in the spin system.
```

### 2.5. **sim\_lcmrawbasis.m**

**USAGE:**

```
[RF,out]=sim_lcmrawbasis(n,sw,Bfield,linewidth,metab,tau1,tau2,addref,makeraw,  
seq)
```

**DESCRIPTION:**

Generate an LCModel .RAW file to be used as an individual metabolite basis spectrum in an LCModel basis set. The relevant characteristics of the acquisition can be specified (pulse sequence, number of points, spectral

width, etc)

**INPUTS:**

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
tau1 = first echo time in [s] (if seq='st', tau1 = TE)  
tau2 = second echo time in [s]. (Used in Press, but not used in SE. (If seq='st', tau2=TM).  
addref = add reference at 0ppm (for use in LCModel makebasis) ['y' or 'n']  
makeraw = make output file for lcmodel ['y' or 'n']  
seq = pulse sequence ['se' for Spin Echo or 'p' for Press]  
metab = one of the following choices  
    'H2O' = Water  
    'Ala' = Alanine  
    'Asp' = Aspartate  
    'PCh' = PhosphoCholine  
    'Cr' = Creatine  
    'PCr' = PhosphoCreatine  
    'GABA' = Gamma-aminobutyric acid (kaiser)  
    'GABA3' = Gamma-aminobutyric acid (de Graaf)  
    'Gln' = Glutamine  
    'Glu' = Glutamate  
    'GSH' = Glutathione  
    'Gly' = Glycine  
    'Ins' = Myo-inositol  
    'Lac' = Lactate  
    'NAA' = N-acetyl aspartate  
    'Scyлло' = Scyllo-inositol  
    'Tau' = Taurine  
    'Asc' = Ascorbate (Vitamin C)  
    'bHB' = beta-Hydroxybutyrate  
    'bHG' = beta-Hydroxyglutarate  
    'Glc' = Glucose  
    'NAAG' = N-acetyl aspartyl glutamate  
    'GPC' = Glycero-phosphocholine  
    'PE' = Phosphoryl ethanolamine  
    'Ser' = Serine

## 2.6. sim\_megapress.m

**USAGE:**

out=sim\_megapress(n,sw,Bfield,linewidth,sys,taus,refoc1Flip,refoc2Flip,editFlip)

**DESCRIPTION:**

Simulate the MEGA-PRESS sequence with instantaneous localization and editing pulses. Provides the ability to specify the flip angle of each refocusing pulse and editing pulse on each spin in the spin system.

**INPUTS:**

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
sys = spin system definition structure  
taus = pulse sequence timing vector:

```

    taus(1)      = time in [ms] from 90 to 1st 180
    taus(2)      = time in [ms] from 1st 180 to 1st edit pulse
    taus(3)      = time in [ms] from 1st edit pulse to 2nd 180
    taus(4)      = time in [ms] from 2nd 180 to 2nd edit pulse
    taus(5)      = time in [ms] from 2nd edit pulse to ADC
    refoc1Flip= array of refoc1 flip angles for each spin in system
    refoc2Flip= array of refoc2 flip angles for each spin in system
    editFlip  = array of editing flip angles for each spin in system

```

## 2.7. **sim\_megapress\_shaped.m**

### USAGE:

```

sim_megapress_shaped(n,sw,Bfield,linewidth,taus,sys,editPulse,editTp,editPh1,e
ditPh2,refPulse,refTp,Gx,Gy,dx,dy,refPh1,refPh2)

```

### DESCRIPTION:

This function simulates the MEGA-PRESS sequence with shaped localization pulses and shaped editing pulses. Enables choice of the timings of all of the rf pulses as well as the choice of the phase of both the editing pulse and the refocusing pulses. This allows phase cycling of the editing and refocusing pulses by repeating simulations with different editing pulse phases, which is necessary to remove phase artefacts from the editing pulses. For the editing pulses, an eight step phase cycling scheme is typically sufficient, where by the first editing pulse is cycled by 0 and 90 degrees, and the second editing pulse is cycled by 0,90,180, and 270 degrees, and all phase cycles should be added together to remove unwanted coherences. For the refocusing pulses, a four step phase cycling scheme is typically sufficient, where both refocusing pulses are phase cycled by 0 and 90 degrees, and the phase are combined in the following way:

```

signal = ([0 90] - [0 0]) + ([90 0] - [90 90]);

```

where, in [X Y], X is the phase of the first refocusing pulse and Y is the phase of the second refocusing pulse

Note that this code only simulates one subspectrum at a time (edit-on or edit-off). The difference spectrum can be obtained by simulating one of each, and then subtracting.

### INPUTS:

```

n          = number of points in fid/spectrum
sw         = desired spectral width in [Hz]
Bfield     = main magnetic field strength in [T]
linewidth  = linewidth in [Hz]
taus(1)    = time in [ms] from 90 to 1st 180
taus(2)    = time in [ms] from 1st 180 to 1st edit pulse
taus(3)    = time in [ms] from 1st edit pulse to 2nd 180
taus(4)    = time in [ms] from 2nd 180 to 2nd edit pulse
taus(5)    = time in [ms] from 2nd edit pulse to ADC
            FOR MEGA-PRESS on SIEMENS SYSTEM:
            taus=[4.545,12.7025,21.7975,12.7025,17.2526];
editPulse  = Editing pulse shape - [Nx3] array with [:,1]=phase,
            [:,2]=amplitude, [:,3]=duration
editTp     = duration of editing pulse in [ms];
sys        = Metabolite spin system definition structure;
editPh1    = the phase of the first editing pulse in [degrees];
editPh2    = the phase of the second editing pulse in [degrees];
refPulse   = Refocusing pulse shape - [Nx3] array with [:,1]=phase,
            [:,2]=amplitude, [:,3]=duration

```



refTp       = duration of refocusing pulse in [ms]  
 Gx         = gradient strength for first selective refocusing pulse [G/cm]  
 Gy         = gradient strength for second selective refocusing pulse [G/cm]  
 dx         = position offset in x-direction (corresponding to first refocusing pulse) [cm]  
 dy         = position offset in y-direction (corresponding to second refocusing pulse) [cm]  
 refPh1      = the phase of the first refocusing pulse in [degrees];  
 refPh2      = the phase of the second refocusing pulse in [degrees];

## 2.8. **sim\_megapress\_shapedEdit.m**

### USAGE:

```
sim_megapress_shapedEdit(n,sw,Bfield,linewidth,taus,sys,editPulse,editTp,editPh1,editPh2)
```

### DESCRIPTION:

This function simulates the MEGA-PRESS sequence with instantaneous localization pulses and shaped editing pulses. Enables choice of the timings of all of the rf pulses as well as the choice of the phase of the editing pulse. This allows phase cycling of the editing pulses by repeating simulations with different editing pulse phases, which is necessary to remove phase artefacts from the editing pulses. For the editing pulses, an eight step phase cycling scheme is typically sufficient, where by the first editing pulse is cycled by 0 and 90 degrees, and the second editing pulse is cycled by 0,90,180, and 270 degrees, and all phase cycles should be added together to remove unwanted coherences.

Note that this code only simulates one subspectrum at a time (edit-on or edit-off). The difference spectrum can be obtained by simulating one of each, and then subtracting.

### INPUTS:

n           = number of points in fid/spectrum  
 sw          = desired spectral width in [Hz]  
 Bfield      = main magnetic field strength in [T]  
 linewidth   = linewidth in [Hz]  
 taus(1)     = time in [ms] from 90 to 1st 180  
 taus(2)     = time in [ms] from 1st 180 to 1st edit pulse  
 taus(3)     = time in [ms] from 1st edit pulse to 2nd 180  
 taus(4)     = time in [ms] from 2nd 180 to 2nd edit pulse  
 taus(5)     = time in [ms] from 2nd edit pulse to ADC  
             FOR MEGA-PRESS on SIEMENS SYSTEM:  
             taus=[4.545,12.7025,21.7975,12.7025,17.2526];  
 editPulse   = Editing pulse shape - [Nx3] array with [:,1]=phase,  
             [:,2]=amplitude, [:,3]=duration  
 editTp      = duration of editing pulse in [ms];  
 sys         = Metabolite spin system definition structure;  
 editPh1     = the phase of the first editing pulse in [degrees];  
 editPh2     = the phase of the second editing pulse in [degrees];

## 2.9. **sim\_megapress\_shapedRefoc.m**

### USAGE:

```
sim_megapress_shapedRefoc(n,sw,Bfield,linewidth,taus,sys,editFlip,refPulse,ref
```

Tp,Gx,Gy,dx,dy,refPh1,refPh2)

DESCRIPTION:

This function simulates the MEGA-PRESS sequence with shaped localization pulses and instantaneous editing pulses. Enables choice of the timings of all of the rf pulses as well as the choice of the phase of both the editing pulse and the refocusing pulses. This allows phase cycling of the editing and refocusing pulses by repeating simulations with different editing pulse phases, which is necessary to remove phase artefacts from the editing pulses. For the editing pulses, an eight step phase cycling scheme is typically sufficient, where by the first editing pulse is cycled by 0 and 90 degrees, and the second editing pulse is cycled by 0,90,180, and 270 degrees, and all phase cycles should be added together to remove unwanted coherences. For the refocusing pulses, a four step phase cycling scheme is typically sufficient, where both refocusing pulses are phase cycled by 0 and 90 degrees, and the phase are combined in the following way:

```
signal = ([0 90] - [0 0]) + ([90 0] - [90 90]);
```

where, in [X Y], X is the phase of the first refocusing pulse and Y is the phase of the second refocusing pulse

Note that this code only simulates one subspectrum at a time (edit-on or edit-off). The difference spectrum can be obtained by simulating one of each, and then subtracting.

INPUTS:

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
taus(1) = time in [ms] from 90 to 1st 180  
taus(2) = time in [ms] from 1st 180 to 1st edit pulse  
taus(3) = time in [ms] from 1st edit pulse to 2nd 180  
taus(4) = time in [ms] from 2nd 180 to 2nd edit pulse  
taus(5) = time in [ms] from 2nd edit pulse to ADC  
FOR MEGA-PRESS on SIEMENS SYSTEM:  
taus=[4.545,12.7025,21.7975,12.7025,17.2526];  
sys = Metabolite spin system definition structure;  
editFlip = vector of editing flip angles in [degrees] at chemical shifts corresponding to 'shifts'.  
refPulse = Refocusing pulse shape - [Nx3] array with [:,1]=phase, [:,2]=amplitude, [:,3]=duration  
refTp = duration of refocusing pulse in [ms]  
Gx = gradient strength for first selective refocusing pulse [G/cm]  
Gy = gradient strength for second selective refocusing pulse [G/cm]  
dx = position offset in x-direction (corresponding to first refocusing pulse) [cm]  
dy = position offset in y-direction (corresponding to second refocusing pulse) [cm]  
refPh1 = the phase of the first refocusing pulse in [degrees];  
refPh2 = the phase of the second refocusing pulse in [degrees];

## 2.10. sim\_onepulse.m

USAGE:

```
out=sim_onepulse(n,sw,Bfield,linewidth,sys)
```

DESCRIPTION:

This function simulates a pulse-acquire experiment with an ideal (instantaneous) excitation pulse and an assumed lorentzian lineshape. The function calls the function 'sim\_Hamiltonian' which produces the free evolution Hamiltonian for the specified number of spins, J and shifts.

**INPUTS:**

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
sys = spin system definition structure

## **2.11. sim\_onepulse\_arbPh.m**

**USAGE:**

out=sim\_onepulse\_arbPh(n,sw,Bfield,linewidth,sys,ph)

**DESCRIPTION:**

This function simulates a pulse-acquire experiment with an ideal (instantaneous) excitation pulse and an assumed lorentzian lineshape. The function calls the function 'sim\_Hamiltonian' which produces the free evolution Hamiltonian for the specified number of spins, J and shifts. This function enables an excitation pulse with an arbitrary phase.

**INPUTS:**

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
sys = spin system definition structure  
ph = excitation pulse phase (in degrees)

## **2.12. sim\_onepulse\_shaped.m**

**USAGE:**

out = sim\_onepulse\_shaped(n,sw,Bfield,linewidth,sys,RF,tp,dwdx,G)

**DESCRIPTION:**

This function simulates the effect of a frequency selective or slice selective excitation, followed immediately by the acquisition window. This is mainly an exercise to see if I can get slice selective excitation working.

Note that when simulating a frequency selective pulse, it is okay to specify only 8 arguments (no gradient needs to be specified). If the 9th argument, G, is specified and is non-zero, then a slice selective pulse is assumed.

**INPUTS:**

n = number of points in fid/spectrum  
sw = desired spectral width in [Hz]  
Bfield = main magnetic field strength in [T]  
linewidth = linewidth in [Hz]  
sys = spin system definition structure  
RF = radiofrequency pulse array [N x 3]. Phase, Amplitude, Duration.  
tp = RF pulse duration in [ms]  
peakB1 = peak B1 amplitude in [kHz]  
dwdx = if simulating a frequency selective pulse, this argument should be

the frequency offset [Hz]. If simulating a slice selective pulse, this argument should be the position offset [cm].  
 G = gradient strength for slice-selective pulse [G/cm];

### 2.13. `sim_press.m`

#### USAGE:

```
out = sim_press(n,sw,Bfield,linewidth,sys,tau1,tau2)
```

#### DESCRIPTION:

This function simulates an ideal PRESS experiment. The function calls the function 'sim\_Hamiltonian.m' which produces the free evolution Hamiltonian for the specified number of spins, J and shifts. This time the individual Iy and Iz are needed as well. This function simulates a spin-echo experiment with echo time tau.

#### INPUTS:

```
n           = number of points in fid/spectrum
sw          = desired spectral width in [Hz]
Bfield      = main magnetic field strength in [T]
linewidth   = linewidth in [Hz]
sys         = spin system definition structure
tau1        = Echo time in [s] of first press Spin Echo
tau2        = Echo time in [s] of second press Spin Echo
```

### 2.14. `sim_press_shaped.m`

#### USAGE:

```
out =
sim_press_shaped(n,sw,Bfield,linewidth,sys,tau1,tau2,RF,tp,dx,dy,Gx,Gy,phCyc1,
phCyc2)
```

#### DESCRIPTION:

This function simulates the PRESS experiment. The excitation is simulated as an instantaneous rotation, and the refocusing pulse is simulated as a shaped rotation.

This code enables the choice of the phase of the refocusing pulses. This enables phase cycling of the refocusing pulses by repeating simulations with different editing pulse phases, which is necessary to remove phase artefacts from the editing pulses. A four step phase cycling scheme is typically sufficient, where both refocusing pulses are phase cycled by 0 and 90 degrees, and the phase are combined in the following way:

```
signal = ([0 90] - [0 0]) + ([90 0] - [90 90]);
```

where, in [X Y], X is the phase of the first refocusing pulse and Y is the phase of the second refocusing pulse.

Finally, this code simulates the spectrum at a given point in space (x,y), given the values of the slice selection gradients (Gx, and Gy). The pulse waveform is assumed to be the same for both refocusing pulses. In order to fully simulate the MEGA-PRESS experiment, you have to run this simulation many times at various points in space (x,y), and then add together the resulting spectra.

#### INPUTS:

n = number of points in fid/spectrum  
 sw = desired spectral width in [Hz]  
 Bfield = main magnetic field strength in [T]  
 linewidth = linewidth in [Hz]  
 sys = spin system definition structure  
 tau1 = echo time 1 in [ms].  
 tau2 = echo time 2 in [ms].  
 RF = radiofrequency pulse array [N x 3]. Phase, Amplitude, Duration.  
 tp = RF pulse duration in [ms]  
 peakB1 = peak B1 amplitude in [kHz]  
 dx = position offset in x-direction (corresponding to first refocusing pulse) [cm]  
 dy = position offset in y-direction (corresponding to second refocusing pulse) [cm]  
 Gx = gradient strength for first selective refocusing pulse [G/cm]  
 Gy = gradient strength for second selective refocusing pulse [G/cm]  
 phCycl1 = initial phase of the first refocusing pulse in [degrees];  
 phCycl2 = initial phase of the second refocusing pulse in [degrees];

## 2.15. sim\_readout.m

### USAGE:

d\_out = sim\_readout(d\_in,H,n,sw,linewidth,rcvPhase,shape)

### DESCRIPTION:

This function simulates an ADC readout of the transverse magnetization during the free evolution of the spin system under the effects of chemical shift and scalar coupling.

### INPUTS:

d\_in = input density matrix structure.  
 H = Hamiltonian operator structure.  
 n = number of readout points  
 sw = spectral width [Hz]  
 linewidth = full width at half maximum of spectral peaks [Hz]  
 rcvPhase = receiver phase [degrees]. Optional. Default = 0 (corresponds to x'-axis readout);  
 shape = line broadening function. Optional,  
         'L' = lorentzian (default)  
         'G' = gaussian

## 2.16. sim\_rotate.m

### USAGE:

d\_out = sim\_rotate(d\_in,H,angle,axis)

### DESCRIPTION:

This function simulates the effect of an ideal (instantaneous) rotation on the density matrix. Used in simulation tools.

### INPUTS:

d\_in = input density matrix structure.  
 H = Hamiltonian operator structure.  
 angle = RF pulse flip angle (degrees). If this value is a scalar, then the same flip angle will be applied to all spins in the system. To apply a different flip angle to the different spins in the system, the angle variable can be a vector of flip angles with

length equal to the H.nspins.  
axis = Axis of rotation ('x', 'y' or 'z'); (A z-rotation technically doesn't correspond to an rf pulse rotation, but it is included here anyway).

## 2.17. sim\_rotate\_arbPh.m

### USAGE:

d\_out = sim\_rotate\_arbPh(d\_in,H,angle,ph)

### DESCRIPTION:

This function simulates the effect of an ideal (instantaneous) rotation on the density matrix. Used in simulation tools. The phase of the rf pulse can be arbitrarily chosen. To achieve an arbitrary phase, this code executes a rotation about z by an angle of -phi (the rf pulse phase), then executes a rotation about x by an angle of "angle" (the flip angle of the pulse), and then executes a rotation back about z by an angle of phi.

### INPUTS:

d\_in = input density matrix structure.  
H = Hamiltonian operator structure.  
angle = RF pulse flip angle (degrees). If this value is a scalar, then the same flip angle will be applied to all spins in the system. To apply a different flip angle to the different spins in the system, the angle variable can be a vector of flip angles with length equal to the H.nspins.  
ph = Phase of rotation (in degrees; ie. 0='x', 90='y');

## 2.18. sim\_runAllMetabs.m

### USAGE:

This script is run simply by editing the input parameters and then clicking "Run".

### DESCRIPTION:

Script to generate simulated basis spectra for all metabolites of interest in the human brain. The script will generate an LCModel format .RAW file for each metabolite basis spectrum, which can then be passed into LCModel's "makebasis" function to generate a complete LCModel basis set.

### INPUTS:

To run this script, edit the following parameters as desired and then click run:

lb = linewidth (Hz)  
np = Spectral points  
sw = Spectral width (Hz)  
Bo = Magnetic Field Strength (Tesla)  
te1 = First PRESS echo time, or SPECIAL echo time (s)  
te2 = Second PRESS echo time (if applicable) (s).  
seq = Pulse sequence ('se'= SPECIAL, 'p'=press, 'st'=steam);  
ref = Add reference peak at 0ppm (used in LCModel, y or n);

## 2.19. sim\_runMegaPressShaped.m

**USAGE:**

This script is run simply by editing the input parameters and then clicking "Run".

**DESCRIPTION:**

This script simulates a MEGA-PRESS experiment with fully shaped editing and refocusing pulses. Phase cycling of both the editing and refocusing pulses is performed. Furthermore, simulations are run at various locations in space to account for the within-voxel spatial variation of the GABA signal. Summation across phase cycles and spatial positions is performed. As a result of the phase cycling and spatially resolved simulations, this code takes a long time to run. Therefore, the MATLAB parallel computing toolbox (parfor loop) was used to accelerate the simulations. Acceleration is currently performed in the direction of the slice selective pulse along the x-direction, but this can be changed. Up to a factor of 12 acceleration can be achieved using this approach. To enable the use of the MATLAB parallel computing toolbox, initialize the multiple worked nodes using "matlabpool size X" where "X" is the number of available processing nodes. If the parallel processing toolbox is not available, then replace the "parfor" loop with a "for" loop.

**INPUTS:**

To run this script, edit the parameters below as desired and then click "run":

refocWaveform	= name of refocusing pulse waveform.
editWaveform	= name of editing pulse waveform.
editOnFreq	= frequency of edit on pulse[ppm]
editOffFreq	= frequency of edit off pulse[ppm]
refTp	= duration of refocusing pulses[ms]
editTp	= duration of editing pulses[ms]
Bfield	= Magnetic field strength in [T]
Npts	= number of spectral points
sw	= spectral width [Hz]
Bfield	= magnetic field strength [Tesla]
lw	= linewidth of the output spectrum [Hz]
thkX	= slice thickness of x refocusing pulse [cm]
thkY	= slice thickness of y refocusing pulse [cm]
x	= vector of X positions to simulate [cm]
y	= vector of y positions to simulate [cm]
taus	= vector of pulse sequence timings [ms]
spinSys	= spin system to simulate
editPhCyc1	= vector of phase cycling steps for 1st editing pulse
[degrees]	
editPhCyc2	= vector of phase cycling steps for 2nd editing pulse
[degrees]	
refPhCyc1	= vector of phase cycling steps for 1st refocusing pulse
[degrees]	
refPhCyc2	= vector of phase cycling steps for 2nd refocusing pulse
[degrees]	

**2.20. sim\_runMegaPressShapedEdit.m****USAGE:**

This script is run simply by editing the input parameters and then clicking "Run".

**DESCRIPTION:**

This script simulates a MEGA-PRESS experiment with fully shaped editing pulses. Phase cycling of editing pulses is performed, and summation across phase cycles is performed.

**INPUTS:**

To run this script, edit the parameters below as desired and then click "run":

```

editWaveform      = name of editing pulse waveform.
editOnFreq        = frequency of edit on pulse[ppm]
editOffFreq       = frequency of edit off pulse[ppm]
editTp           = duration of editing pulses[ms]
Npts             = number of spectral points
sw              = spectral width [Hz]
Bfield           = magnetic field strength [Tesla]
lw              = linewidth of the output spectrum [Hz]
taus            = vector of pulse sequence timings [ms]
spinSys          = spin system to simulate
editPhCyc1       = vector of phase cycling steps for 1st editing pulse
                  [degrees]
editPhCyc2       = vector of phase cycling steps for 2nd editing pulse
                  [degrees]

```

**2.21. sim\_runMegaPressShapedRefoc.m****USAGE:**

This script is run simply by editing the input parameters and then clicking "Run".

**DESCRIPTION:**

This script simulates a MEGA-PRESS experiment with fully shaped refocusing pulses. Phase cycling of the refocusing pulses is performed and simulations are run at various locations in space to account for the within-voxel spatial variation of the GABA signal. Summation across spatial positions is performed. As a result of the phase cycling and spatially resolved simulations, this code takes a long time to run. Therefore, the MATLAB parallel computing toolbox (parfor loop) is used to accelerate the simulations. Acceleration is currently performed in the direction of the slice selective pulse along the x-direction, but this can be changed. Up to a factor of 12 acceleration can be achieved using this approach. To enable the use of the MATLAB parallel computing toolbox, initialize the multiple worked nodes using "matlabpool size X" where "X" is the number of available processing nodes. If the parallel processing toolbox is not available, then replace the "parfor" loop with a "for" loop.

**INPUTS:**

To run this script, edit the parameters below as desired and then click "run":

```

refocWaveform     = name of refocusing pulse waveform.
refTp            = duration of refocusing pulses[ms]
Npts             = number of spectral points
sw              = spectral width [Hz]
Bfield           = magnetic field strength [Tesla]
lw              = linewidth of the output spectrum [Hz]
thkX             = slice thickness of x refocusing pulse [cm]
thkY             = slice thickness of y refocusing pulse [cm]
x               = vector of X positions to simulate [cm]
y               = vector of y positions to simulate [cm]
taus            = vector of pulse sequence timings [ms]
spinSys          = spin system to simulate
editFlipON       = vector of edit-ON pulse flip angles for each spin in
                  spin system.
editFlipOFF      = vector of edit-OFF pulse flip angles for each spin in
                  spin system.
refPhCyc1       = vector of phase cycling steps for 1st refocusing
                  pulse [degrees]

```



refPhCyc2                      = vector of phase cycling steps for 2nd refocusing pulse [degrees]

## 2.22. sim\_shapedRF.m

### USAGE:

d\_out = sim\_shapedRF(d\_in,H,RFstruct,flipAngle,phase,grad,pos)

### DESCRIPTION:

This function simulates the effect of a shaped rf pulse on the density matrix. The temporal shape of the refocussing pulses is modelled as a series of N instantaneous rotations about the effective RF field, where N is the number of time points in the RF waveform. The instantaneous effective RF field can be an arbitrary vector, and can be represented in polar coordinates as  $B(\text{Beff}, \alpha, \text{zeta})$ , where Beff is the magnitude field, alpha is the polar angle (the angle between the transverse plane and the effective B field), and zeta is the azimuthal angle, which is given by the phase of the RF). Rotation about the effective B-field is achieved by a composite rotation: Rotate about Y by  $-\alpha$ , rotate about Z by  $-\text{zeta}$ , then rotate about X by  $2\pi \cdot \gamma \cdot \text{Beff} \cdot \text{dt}$ , then rotate back about Z by zeta and back about Y by alpha.

### INPUTS:

d\_in            = input density matrix structure.  
H              = Hamiltonian operator structure.  
RF             = Radiofrequency pulse. This can be the filename of a Siemens .pta file, or an RF pulse definition structure (obtained using rf\_init.m)  
Tp             = Pulse duration in [ms];  
flipAngle      = RF pulse flip angle [degrees].  
phase          = Phase of RF pulse [degrees]. Optional. Default = 0 (x'-axis. 90 degrees corresponds to +y' axis)  
grad           = Gradient strength [G/cm]. Optional (for slice selective pulses only).  
pos            = Position of spins relative to voxel centre [cm]. Optional (for slice selective pulses only).

## 2.23. sim\_spinecho.m

### USAGE:

out = sim\_spinecho(n,sw,Bfield,linewidth,sys,tau)

### DESCRIPTION:

This function simulates a spin-echo experiment with instantaneous RF pulses.

### INPUTS:

n              = number of points in fid/spectrum  
sw             = desired spectral width in [Hz]  
Bfield        = main magnetic field strength in [T]  
linewidth     = linewidth in [Hz]  
sys            = spin system definition structure  
tau            = echo time in [s]

## 2.24. sim\_spinecho\_shaped.m

**USAGE:**

```
out = sim_spinecho_shaped(n,sw,Bfield,linewidth,sys,tau,RF,Tp,grad,pos,ph)
```

**DESCRIPTION:**

This function simulates a spin-echo experiment with instantaneous RF pulses.

**INPUTS:**

```
n          = number of points in fid/spectrum
sw         = desired spectral width in [Hz]
Bfield     = main magnetic field strength in [T]
linewidth  = linewidth in [Hz]
sys        = spin system definition structure
tau        = echo time in [s]
```

**2.25. sim\_spoil.m****USAGE:**

```
d_out = sim_spoil(d_in,H,angle)
```

**DESCRIPTION:**

This function simulates the effect of a rotation about the z-axis.

**INPUTS:**

```
d_in      = input density matrix structure.
H         = Hamiltonian operator structure.
angle     = Spoil angle (degrees).
```

**2.26. sim\_steam.m****USAGE:**

```
out = sim_steam(n,sw,Bfield,linewidth,sys,te,tm)
```

**DESCRIPTION:**

Simulate the STEAM sequence using ideal (instantaneous) RF pulses. To remove unwanted coherences, a 4 step phase cycle is automatically performed, with the first and third rf pulses being cycled by 0, 90 180, and 270 degrees. THIS CODE IS NOT TESTED. RESULTS MAY NOT BE ACCURATE!!

**INPUTS:**

```
n          = number of points in fid/spectrum
sw         = desired spectral width in [Hz]
Bfield     = main magnetic field strength in [T]
linewidth  = linewidth in [Hz]
sys        = spin system definition structure
te         = echo time in [s]
tm         = mixing time in [s]
```

**2.27. sim\_steam\_gradSim.m****USAGE:**

Script can be run by pressing "run".

**DESCRIPTION:**

This function runs the `sim_steam_spoil` function multiple times with different spoiler gradient intensities. The result is a spoiled STEAM sequence.

**INPUTS:**

Initialize the following variables and then click "run":

`spinsys` = Spin system.  
`TE` = Echo time [s].  
`TM` = Mixing time [s].  
`N` = Number of 'phase cycles'

## **2.28. writelcmraw.m**

**USAGE:**

`RF=writelcmraw(data_struct,outfile,metab);`

**DESCRIPTION:**

Take a simulated metabolite basis spectrum in matlab structure format, and output it into LCModel RAW format to be used in an LCModel basis spectrum.

**INPUTS:**

`data_struct` = simulated metabolite basis spectrum in matlab structure format.  
`outfile` = name of the output .RAW file.  
`metab` = Abbreviated name of the metabolite (ie. 'Cr', 'Glu', etc.)

## **3. RF Pulse Tools**

### **3.1. rf\_blochSim.m**

**USAGE:**

`[mv,sc]=rf_blochSim(RF,tp,peakB1,sc);`

**DESCRIPTION:**

Perform a bloch simulation of an RF pulse.

**INPUTS:**

`RF` = RF pulse definition structure  
`tp` = pulse duration in [ms]  
`sc` = Frequency span in [kHz] (optional)  
`peakB1` = Peak B1 amplitude in [kHz] (optional)

### **3.2. rf\_dualBand.m**

**USAGE:**

`[rf,AMPINT]=rf_dualBand(tp,df,n,bw,ph,shft)`

**DESCRIPTION:**

Creates an n-point dual banded gaussian inversion RF pulse with duration

tp(ms). The first band will be at  $f=0\text{Hz}$  and the second band will be at  $df\text{ Hz}$ .  
Bw is the bandwidth of the two selection bands in Hz.

**INPUTS:**

tp           = pulse duration in ms.  
df           = frequency of 2nd gaussian band [Hz].  
n            = number of points in rf waveform.  
bw           = bandwidth of both selection bands [Hz].  
ph           = phase of the second gaussian.  
shift        = frequency shift applied to both bands.

### **3.3. rf\_freqshift.m**

**USAGE:**

RF\_shift=rf\_freqshift(RF,Tp,f);

**DESCRIPTION:**

Apply a frequency shift to an RF pulse.

**INPUTS:**

RF           = RF pulse definition structure.  
Tp           = duration of the rf pulse in [ms].  
f            = amount that you would like to frequency shift the rf pulse in [Hz].

### **3.4. rf\_gauss.m**

**USAGE:**

[rf,AMPINT]=rf\_gauss(tp,df,n,bw);

**DESCRIPTION:**

Create an n point gaussian rf waveform. Waveform can be converted in to siemens pta file using rf\_writepta.m or into varian/agilent rf file using rf\_writerf.m.

**INPUTS:**

tp           = duration of rf pulse in ms.  
df           = frequency of Gaussian pulse in Hz. 0 frequency will correspond to the reference frequency of the rf transmitter.  
n            = number of points in the rf waveform.  
bw           = FWHM of the Gaussian inversion profile in the frequency domain (Hz).

### **3.5. rf\_hs.m**

**USAGE:**

[rf,FM,scan,mvector]=rf\_hs(outfile,N,n,tbw,Tp,trunc,thk)

**DESCRIPTION:**

This function creates any desired HS pulse. N is the number of steps, n is the order of the HS pulse, tbw is the time bandwidth product of the pulse, Tp is the duration of the pulse and thk is the desired thickness of the pulse.

INPUTS:

outfile	= name of output rf file.
N	= Number of points in RF waveform.
n	= order of the HS pulse.
tbw	= Time bandwidth product.
TP	= Duration of the RF pulse (ms).
trunc	= Truncation of the amplitude modulation function.
thk	= thickness of the slice selective pulse.

### 3.6. rf\_loadwaveform.m

USAGE:  
[RF\_struct]=rf\_loadwaveform(filename,type,f0);

DESCRIPTION:  
Initialize an RF pulse structure to contain an RF Pulse waveform as well as its accompanying header information. This function finds the time-bandwidth product (tbw) and the time-w1 product (tw1) of the pulse, and stores this information in the header fields of the output RF structure.

INPUTS:

filename	= filename of RF pulse waveform text file. Can be in Siemens format (.pta) or Varian/Agilent format (.RF).
type	= Excitation ('exc'), Refocusing ('ref') or Inversion ('inv')
f0	= centre frequency of the rf pulse [Hz]. Optional. Default=0.

### 3.7. rf\_readRF.m

USAGE:  
[rf,info]=rf\_readRF(filename)

DESCRIPTION:  
Read a Varian/Agilent .RF file into matlab. The resulting RF matrix will have 3 columns specifying magnitude, phase, and duration. This function simply calls Martyn Klassen's readrfvnmr.m function.

INPUTS:

filename	= filename of the .RF file to read in.
----------	--

### 3.8. rf\_readpta.m

USAGE:  
[rf,info]=rf\_readpta(filename)

DESCRIPTION:  
Read a Siemens .pta file into matlab. The resulting RF matrix will have 2 columns specifying magnitude and phase.

INPUTS:

filename	= filename of the .pta file to read in.
----------	---

### **3.9. rf\_resample.m**

**USAGE:**

RF\_out=rf\_resample(RF\_in,N);

**DESCRIPTION:**

Resample the input RF pulse into a new waveform with N discrete points.

**INPUTS:**

RF           = RF pulse definition structure  
N            = Number of points in new RF waveform

### **3.10. rf\_writepta.m**

**USAGE:**

RF=rf\_writepta(rf,outfile);

**DESCRIPTION:**

Write a matlab RF pulse structure (containing 3 x N waveform array field with rf.waveform(1,:)= phase, rf.waveform(2,:)=amplitude, and rf.waveform(3,:)=timestep), to a siemens format .pta file.

**INPUTS:**

rf           = matlab RF pulse.  
outfile      = name of the output .pta file to be written.

### **3.11. rf\_writerf.m**

**USAGE:**

RF=rf\_writerf(rf,outfile);

**DESCRIPTION:**

Write a matlab RF pulse structure (containing 3 x N waveform array field with rf.waveform(1,:)= phase, rf.waveform(2,:)=amplitude, and rf.waveform(3,:)=timestep), to a varian/agilent format .RF file.

**INPUTS:**

rf           = matlab RF pulse.  
outfile      = name of the output .RF file to be written.

## **4. Processing Tools**

### **4.1. addphase.m**

**USAGE:**

PhasedSpecs=addphase(specs,AddedPhase);

**DESCRIPTION:**

Add equal amount of complex phase to each point of a vector.

INPUTS:  
specs = Input vector.  
AddedPhase = Amount of phase (degrees) to add.

#### 4.2. addphase1.m

USAGE:  
PhasedSpecs=addphase1(specs,ppm,timeShift,ppm0,B0);

DESCRIPTION:  
Add first order phase to a spectrum (added phase is linearly dependent on frequency).

INPUTS:  
specs = input vector  
ppm = frequency scale (ppm) corresponding to the specs vector  
timeShift = This defines the amount of 1st order phase shift by specifying the equivalent horizontal shift (in seconds) in the time domain.  
ppm0 = The frequency "origin" (ppm) of the 1st order phase shift. (this frequency will undergo 0 phase shift).  
B0 = The main magnetic field strength (needed since ppm depends on B0)

#### 4.3. getLW.m

USAGE:  
[FWHM]=getLW(in,zpfactor,Refppmmin,Refppmmax);

DESCRIPTION:  
Estimates the linewidth of a reference peak in the spectrum. Two methods are used to estimate the linewidth: 1. FWHM is measured by simply taking the full width at half max of the reference peak. 2. The FWHM is measured by fitting the reference peak to a lorentzian lineshape and determine the FWHM of the best fit. The output FWHM is given by the average of these two measures.

INPUTS:  
in = input spectrum in structure format.  
zpfactor = zero-padding factor (used for method 1.)  
Refppmmin = Min of frequency range (ppm) in which to search for reference peak.  
Refppmmax = Max of frequency range to (ppm) in which search for reference peak

#### 4.4. getLWandSNR.m

USAGE:  
[FWHMmean,SNRmean] = getLWandSNR(in);

DESCRIPTION:  
Calculate the linewidth and SNR of a spectrum.

**INPUTS:**

in = input data in matlab structure format

#### **4.5. getSNR.m**

**USAGE:**

```
[SNR]=getSNR(in,NAAppmmin,NAAppmmax,noiseppmmin,noiseppmmax);
```

**DESCRIPTION:**

Find the SNR of the NAA peak in a spectrum.

**INPUTS:**

in = input data in matlab structure format  
NAAppmmin = min of frequency range in which to search for NAA peak.  
NAAppmmax = max of frequency range in which to search for NAA peak.  
noiseppmmin = min of frequency range in which to measure noise.  
noiseppmmax = max of frequency range in which to measure noise.

#### **4.6. op\_ISIScombine.m**

**USAGE:**

```
out=op_ISIScombine(in,addInd);
```

**DESCRIPTION:**

Combine dimensions corresponding to ISIS on/off acquisitions to produce fully localized MRS volumes. Mostly used for MEGA-SPECIAL data to combine the ISIS subspectra but not the Edit-on/edit-off subspectra.

**INPUTS:**

in = input data in matlab structure format.  
addInd = (optional) If add==1, then row indices [1 2] and [3 4] will be added. Otherwise, row indices [1 2] and [3 4] will be subtracted.

#### **4.7. op\_addNoise.m**

**USAGE:**

```
[out,noise]=op_addNoise(in,sdnoise,noise);
```

**DESCRIPTION:**

Add noise to a spectrum (useful for generating simulated data). Normally distributed random noise is added to both the real and imaginary parts of the data. Real and imaginary noise parts are uncorrelated.

**INPUTS:**

in = Input data in matlab structure format.  
sdnoise = Standard deviation of the random noise to be added in the time domain.  
noise = (optional) Specific noise kernel to be added (if specified, sdnoise variable is ignored).



#### 4.8. **op\_addScans.m**

**USAGE:**

```
out=op_addScans(in1,in2,subtract);
```

**DESCRIPTION:**

Add or subtract two scans together.

**INPUTS:**

```
in1          = First spectrum to add (in matlab structure format)
in2          = Second spectrum to add (also in matlab structure format).
subtract     = (optional). Add or subtract? (0 = add, 1=subtract). Default=0;
```

#### 4.9. **op\_addphase.m**

**USAGE:**

```
out=op_addphase(in,ph0,ph1,ppm0,suppressPlot);
```

**DESCRIPTION:**

add zero and first order phase to a spectrum.

**INPUTS:**

```
in           = input spectrum in matlab structure format
ph0          = zero order phase to add (degrees)
ph1          = 1st order phase to add (in seconds);
ppm0         = (optional) frequency reference point. Default = 4.65;
suppressPlot = (optional) Boolean to suppress plots. Default = 0;
```

#### 4.10. **op\_addphaseSubspec.m**

**USAGE:**

```
out=op_addphaseSubspec(in,ph);
```

**DESCRIPTION:**

Add zero order phase to one of the subspectra in a dataset. For example, the edit-on spectrum of a mega-press acquisition.

**INPUTS:**

```
in   = Input spectrum in matlab structure format.
ph   = Phase (in degrees) to add to the second subspectrum.
```

#### 4.11. **op\_addrcvrs.m**

**USAGE:**

```
[out,fids_presum,specs_presum,ph,sig]=op_addrcvrs(in,point,mode,coilcombos);
```

**DESCRIPTION:**

Perform weighted coil recombination for MRS data acquired with a receiver coil array.

**INPUTS:**

```
in           = input spectrum in matlab structure format.
point        = point of fid to use for phase estimation.
mode         = -'w' performs amplitude weighting of channels based on the
              maximum signal of each coil channel.
```

- 'h' performs amplitude weighting of channels based on the maximum signal of each coil channel divided by the square of the noise in each coil channel (as described by Hall et al. Neuroimage 2014).

coilcombos = (optional) The predetermined coil phases and amplitudes as generated by the op\_getcoilcombos.m function. If this argument is provided, the 'point', and 'mode', arguments will be ignored.

#### 4.12. op\_alignAverages.m

**USAGE:**

```
[out,fs,phs]=op_alignAverages(in,tmax,avg,initPars);
```

**DESCRIPTION:**

Perform spectral registration in the time domain to correct frequency and phase drifts. As described in Near et al. Frequency and phase drift correction of magnetic resonance spectroscopy data by spectral registration in the time domain. Magn Reson Med. 2014 Jan 16. doi: 10.1002/mrm.25094. [Epub ahead of print]

**INPUTS:**

in = Input data structure.  
tmax = Maximum time (s) in time domain to use for alignment.  
avg = Align averages to the average of the averages? (y or n)  
initPars = (Optional) Initial fit parameters [freq(Hz), phase(degrees)].  
Default=[0,0];

#### 4.13. op\_alignAverages\_fd.m

**USAGE:**

```
[out,fs,phs]=op_alignAverages_fd(in,minppm,maxppm,tmax,avg,initPars);
```

**DESCRIPTION:**

Perform time-domain spectral registration using a limited range of frequencies to correct frequency and phase drifts. As described in Near et al. Frequency and phase drift correction of magnetic resonance spectroscopy data by spectral registration in the time domain. Magn Reson Med. 2014 Jan 16. doi: 10.1002/mrm.25094. [Epub ahead of print]

**INPUTS:**

in = Input data structure.  
Minppm = Minimum of frequency range (ppm).  
Maxppm = Maximum of frequency range (ppm).  
tmax = Maximum time (s) in time domain to use for alignment.  
avg = Align averages to the average of the averages? (y or n)  
initPars = (Optional) Initial fit parameters [freq(Hz), phase(degrees)].  
Default=[0,0];

#### 4.14. op\_alignScans.m

**USAGE:**

```
[out,ph,frq]=op_alignScans(in,in1,tmax,mode);
```

**DESCRIPTION:**

Use spectral registration to align two separate scans (align in to in1);

**INPUTS:**

in = input (spectrum to be registered)  
in1 = base (spectrum that the input is to be registered to).  
tmax = Maximum time (s) in time domain to use for registration.  
mode = (optional) 'f' - Frequency align only  
          'p' - Phase align only  
          'fp or pf' - Frequency and phase align (default)

#### **4.15. op\_alignScans\_fd.m**

**USAGE:**

[out,ph,frq]=op\_alignScans\_fd(in,in1,fmin,fmax,tmax,mode);

**DESCRIPTION:**

Use spectral registration on a limited frequency range to align two separate MRS datasets (align in to in1);

**INPUTS:**

in = input (spectrum to be registered to the base)  
in1 = base (spectrum that the input is to be registered to).  
fmin = Minimum frequency for spectral alignment (ppm).  
fmax = Maximum frequency for spectral alignment (ppm).  
tmax = Maximum time (s) in time domain to use for registration.  
mode = (optional) 'f' - Frequency align only  
          'p' - Phase align only  
          'fp or pf' - Frequency and phase align (default)

#### **4.16. op\_alignrcvrs.m**

**USAGE:**

[out,ph,sig]=op\_alignrcvrs(in,point,mode,coilcombos);

**DESCRIPTION:**

phase align the receiver channels without combining them.

**INPUTS:**

in = input spectrum in matlab structure format.  
point = Index of point in time domain to use for phase reference.  
mode = -'w' performs amplitude weighting of channels based on the maximum signal of each coil channel.  
          -'h' performs amplitude weighting of channels based on the maximum signal of each coil channel divided by the square of the noise in each coil channel (as described by Hall et al. Neuroimage 2014).  
Coilcombos = (optional) The predetermined coil phases and amplitudes as generated by the op\_getcoilcombos.m function. If this argument is provided, the 'point', and 'mode', arguments will be ignored.

#### **4.17. op\_ampScale.m**

USAGE:  
out=op\_ampScale(in,A);

DESCRIPTION:  
Scale the amplitude of a spectrum by factor A.

INPUTS:  
in = input data in matlab structure format  
A = Amplitude scaling factor.

#### **4.18. op\_averaging.m**

USAGE:  
out=op\_averaging(in);

DESCRIPTION:  
Combine the averages in a scan by adding the averages together. For historical reasons, this function does not divide by the number of averages, and therefore it is not really taking the average of the averages. Dividing by the number of averages is typically performed in a separate step using op\_avgNormalize.

INPUTS:  
in = input data in matlab structure format.

#### **4.19. op\_avgNormalize.m**

USAGE:  
out=op\_avgNormalize(in,scaleFactor);

DESCRIPTION:  
Divide by the number of averages. This is typically performed after op\_averaging.m to obtain a truly "averaged spectrum", since op\_averaging.m adds the averages together but does not divide by the number of averages.

INPUTS:  
in = input data in matlab structure format.  
scaleFactor = (optional) Factor to divide by. default = in.averages.

#### **4.20. op\_combinesubspecs.m**

USAGE:  
out=op\_combinesubspecs(in,mode);

DESCRIPTION:  
Combine the subspectra in an acquisition either by addition or subtraction.

INPUTS:  
in = input data in matlab structure format.  
mode = "-diff" adds the subspectra together (this is counter intuitive, but the reason is that many "difference editing" sequences use

phase cycling of the readout ADC to achieve "subtraction by addition".  
-"summ" performs a subtraction of the subspectra.

#### **4.21. op\_complexConj.m**

**USAGE:**

out=op\_complexConj(in)

**DESCRIPTION:**

take the complex conjugate of the data;

**INPUTS:**

in = Input data in matlab structure format.

#### **4.22. op\_concatAverages.m**

**USAGE:**

out=op\_concatAverages(in1,in2);

**DESCRIPTION:**

Concatenate two scans along the averages dimension. Two scans with 50 averages each will now look like a single scan with 100 averages.

**INPUTS:**

in1 = first input in matlab structure format.

in2 = second input in matlab structure format.

#### **4.23. op\_creFit.m**

**USAGE:**

parsFit=op\_creFit(in,ph0,ph1);

**DESCRIPTION:**

Perform a Lorentzian lineshape fit to the creatine resonance in a brain proton MRS dataset.

**INPUTS:**

in = input data in matlab stucture format.

ph0 = zero order phase to add to the input data.

ph1 = 1st order phase to add to the input data.

#### **4.24. op\_dccorr.m**

**USAGE:**

out=op\_dccorr(in,mode,var);

**DESCRIPTION:**

Do a DC Correction on the data. This method is a frequency domain operation.

**INPUTS:**

In = input data in matlab structure format.

mode = Point('p') or Value('v'). In point mode, the DC offset is calculated automatically at a specific point in the spectrum. In value mode, the user has to provide the value of the desired DC offset.

Var = If mode is 'p', then 'var' is the index of the spectral point that you wish to use to calculate the DC offset. If mode is 'v', then 'var' is the value of the dc offset correction that you wish to employ.

#### **4.25. op\_downsamp.m**

USAGE:

```
out=op_downsamp(in,dsFactor);
```

DESCRIPTION: Change the time domain sampling rate of a spectrum by a factor of 'dsFactor'. Nearest neighbour interpolation is performed by default.

INPUTS:

in = input data in matlab structure format.

dsFactor = factor by which to divide the sampling rate of the fid.

#### **4.26. op\_ecc.m**

USAGE:

```
[out,outw]=op_ecc(in,inw);
```

DESCRIPTION:

Perform an eddy current correction by estimating any non-linearity in the phase of the water unsuppressed data in the time domain and applying the appropriate correction to both the water suppressed and water unsuppressed data.

INPUTS:

in = water suppressed input data in matlab structure format.

in2 = water unsuppressed input data in matlab structure format.

#### **4.27. op\_fddccorr.m**

USAGE:

```
out=op_fddccorr(in,npts);
```

DESCRIPTION:

Correct and DC offset in the frequency domain. This is equivalent to a vertical shift in the frequency domain. The required vertical shift is calculated by taking the average of the first and last "NPTS" points in the frequency domain. This requires that those points are in the noise floor.

INPUTS:

in = input data in matlab structure format.

npts = number of points at both edges of the frequency domain that will be used for estimation of the DC offset of the spectrum.

#### **4.28. op\_filter.m**

**USAGE:**

```
out=op_filter(in,lb);
```

**DESCRIPTION:**

Perform line broadening by multiplying the time domain signal by an exponential decay function.

**INPUTS:**

```
in      = input data in matlab structure format.  
lb      = line broadening factor in Hz.
```

#### **4.29. op\_freqAlignAverages.m**

**USAGE:**

```
[out,fs]=op_freqAlignAverages(in,tmax,avg,initPars);
```

**DESCRIPTION:**

Perform spectral registration in the time domain using only frequency adjustment (no phase adjustment).

**INPUTS:**

```
in          = Input data structure.  
tmax        = Maximum time (s) in time domain to use for alignment.  
avg         = Align averages to the average of the averages? (y or n)  
initPars    = (Optional) Initial fit parameters [freq(Hz), phase(degrees)].  
              Default=[0,0];
```

#### **4.30. op\_freqAlignAverages\_fd.m**

**USAGE:**

```
[out,fs,phs]=op_freqAlignAverages_fd(in,minppm,maxppm,tmax,avg,initPars);
```

**DESCRIPTION:**

Perform time-domain spectral registration using a limited range of frequencies and using only frequency adjustment (no phase adjustment).

**INPUTS:**

```
in          = Input data structure.  
minppm      = Minimum of frequency range (ppm).  
maxppm      = Maximum of frequency range (ppm).  
tmax        = Maximum time (s) in time domain to use for alignment.  
avg         = Align averages to the average of the averages? (y or n)  
initPars    = (Optional) Initial fit parameters [freq(Hz), phase(degrees)].  
              Default=[0,0];
```

#### **4.31. op\_freqrange.m**

**USAGE:**

```
out=op_freqrange(in,ppmmin,ppmmax);
```

**DESCRIPTION:**

Output only a specified frequency range of the input spectrum.

**INPUTS:**

in           = input data in matlab structure format.  
ppmmin       = minimum extent of frequency range in ppm.  
ppmmax       = maximum extent of frequency range in ppm.

#### **4.32. op\_freqshift.m**

**USAGE:**

out=op\_freqshift(in,f);

**DESCRIPTION:**

Apply a frequency shift to the input spectrum by 'f' Hz.

**INPUTS:**

in           = input data in matlab structure format.  
f            = frequency shift to apply (in Hz).

#### **4.33. op\_freqshiftSubspec.m**

**USAGE:**

out=op\_freqshiftSubspec(in,f);

**DESCRIPTION:**

Apply a frequency shift to only one of the sub-spectra in a dataset. This is used to minimize subtraction artefacts from MEGA\_PRESS data, for instance.

**INPUTS:**

in           = input data in matlab structure format.  
f            = frequency shift to apply to subspectrum (in Hz).

#### **4.34. op\_gaussianPeak.m**

**USAGE:**

out=op\_gaussianPeak(n,sw,Bo,lw,ppm0,amp);

**DESCRIPTION:**

Generate a noiseless spectrum containing a single gaussian peak with desired parameters (frequency, amplitude, linewidth, etc.).

**INPUTS:**

n           = Number of points in spectrum.  
sw          = spectral width of spectrum (Hz).  
Bo          = Magnetic field strength (Tesla).  
lw          = Linewidth of gaussian peak (Hz).  
ppm0        = Frequency of gaussian peak (ppm).  
amp         = Amplitude of gaussian peak.

#### **4.35. op\_getcoilcombos.m**

**USAGE:**

coilcombos=op\_getcoilcombos(file\_or\_struct,point);



**DESCRIPTION:**

This function finds the relative coil phases and amplitudes. Coil phases are found by determining the phase and amplitude of the pointth point in the time domain.

**INPUTS:**

file\_or\_struct = this function will accept either a string filename or the name of a structure. If the input is a string, the program will read in the data corresponding to that filename. If the input is a structure, it will operate on that structure.

point = The index of the datapoint in the fid that is used for determination of signal intensity and phase.

#### **4.36. op\_getcoilcombos\_specReg.m**

**USAGE:**

```
coilcombos=op_getcoilcombos_specReg(file_or_struct,tmin,tmax,point);
```

**DESCRIPTION:**

This function finds the relative coil phases and amplitudes. Coil phases are found by fitting in the time domain. (In the original op\_getcoilcombos, the phases were determined simply by observation of the pointth point in the time domain). The "Base" receiver channel is determined as the one with the highest signal (all other coil channels will be registered to that channel).

**INPUTS:**

file\_or\_struct = this function will accept either a string filename or the name of a structure. If the input is a string, the program will read in the data corresponding to that filename. If the input is a structure, it will operate on that structure.

tmin = The earliest timepoint in the fid to be used for spectral registration.

tmax = The latest timepoint in the fid to be used for spectral registration.

point = The index of the datapoint in the fid that is used for determination of signal intensity.

#### **4.37. op\_leftshift.m**

**USAGE:**

```
out=op_leftshift(in,ls);
```

**DESCRIPTION:**

Remove leading datapoints from the fid to get rid of 1st order phase errors.

**INPUTS:**

in = input data in matlab structure format.

ls = number of points to remove from the beginning of the fid.

#### **4.38. op\_loadjmrui.m**

USAGE:  
out=op\_loadjmrui(filename);

DESCRIPTION:  
Load a jMRUI text file into matlab structure format.

INPUTS:  
filename = filename of the jMRUI txt file.

#### **4.39. op\_loadlcmdetail.m**

USAGE:  
[metabs,corrMatrix]=op\_loadlcmdetail(filename);

DESCRIPTION:  
This function loads in the "detailed output" of LCModel and returns the matrix of metabolite correlation coefficients.

INPUTS:  
filename = Filename of the lcmodel detailed output file.

#### **4.40. op\_loadspec\_IMA.m**

USAGE:  
out=op\_loadspec\_IMA(filename,Bo,spectralwidth);

DESCRIPTION:  
Loads a siemens .IMA file into matlab structure format.

INPUTS:  
filename = Filename of Siemens .IMA file to load.  
Bo = Field strength (Tesla).  
spectralwidth = spectral width of the input spectrum (Hz).

#### **4.41. op\_loadspec\_twix.m**

USAGE:  
out=op\_loadspec\_twix(filename);

DESCRIPTION:  
Reads in siemens twix raw data (.dat file) using the mapVBVD.m and twix\_map\_obj.m functions from Philipp Ehses (philipp.ehses@tuebingen.mpg.de).

op\_loadspec\_twix outputs the data in structure format, with fields corresponding to time scale, fids, frequency scale, spectra, and header fields containing information about the acquisition. The resulting matlab structure can be operated on by the other functions in this MRS toolbox.

INPUTS:  
filename = filename of Siemens twix data to load.

#### **4.42. op\_loadspec\_varian.m**

**USAGE:**

```
out=op_loadspec_varian(filename);
```

**DESCRIPTION:**

Reads in varian .fid data using the readfid.m and readprocpa.m functions from Martyn Klassen (mklassen@robarts.ca).

op\_loadspec\_varian outputs the data in structure format, with fields corresponding to time scale, fids, frequency scale, spectra, and header fields containing information about the acquisition. The resulting matlab structure can be operated on by the other functions in this MRS toolbox.

**INPUTS:**

filename = filename of Varian .fid data to load.

#### **4.43. op\_lorentz.m**

**USAGE:**

```
y=op_lorentz(pars,ppm);
```

**DESCRIPTION:**

Generate a parametrized lorentzian peak. This function is fed into fitting tools to enable fitting of peaks using lorentzian lineshapes.

**INPUTS:**

pars = The parameters of the Lorentzian function. This is a five element vector consisting of the following fields:

- Amplitude,
- FWHM, (In Hz)
- Centre Freq, (In ppm)
- baseline offset, (in Amplitude units)
- Phase shift; (in degrees)

ppm = frequency axis vector (in ppm);

#### **4.44. op\_lorentz\_linbas.m**

**USAGE:**

```
y=op_lorentz_linbas(pars,ppm);
```

**DESCRIPTION:**

Generate a parametrized lorentzian peak with a linearly sloping baseline. This function is fed into fitting tools to enable fitting of peaks using lorentzian lineshapes.

**INPUTS:**

pars = The parameters of the Lorentzian function. This is a six element vector consisting of the following fields:

- [ Amplitude,
- FWHM, (In Hz)
- Centre Freq, (In ppm)
- baseline slope, (in Amplitude units per ppm)
- baseline offset, (in Amplitude units)
- Phase shift]; (in degrees)

ppm = frequency axis vector (in ppm);

#### **4.45. op\_lorentzianPeak.m**

**USAGE:**

```
out=op_lorentzianPeak(n,sw,Bo,lw,ppm0,amp);
```

**DESCRIPTION:**

Generate a noiseless spectrum containing a single lorentzian peak with desired parameters (frequency, amplitude, linewidth, etc.).

**INPUTS:**

n        = Number of points in spectrum.  
sw       = spectral width of spectrum (Hz).  
Bo       = Magnetic field strength (Tesla).  
lw       = Linewidth of gaussian peak (Hz).  
ppm0     = Frequency of gaussian peak (ppm).  
amp      = Amplitude of gaussian peak.

#### **4.46. op\_makeFreqDrift.m**

**USAGE:**

```
[out,fDrift]=op_makeFreqDrift(in,totalDrift,noise);
```

**DESCRIPTION:**

Add frequency drift to a dataset containing multiple averages. This is generally used to generate simulated datasets with phase drift.

**INPUTS:**

in                = input data in matlab structure format.  
totalDrift       = total amount of frequency drift (in Hz) to add over the whole scan.  
noise            = the standard deviation of noise to add to the frequency drift function.

#### **4.47. op\_makePhaseDrift.m**

**USAGE:**

```
[out,phDrift]=op_makePhaseDrift(in,totalDrift,noise);
```

**DESCRIPTION:**

Add phase drift to a dataset containing multiple averages. This is generally used to generate simulated datasets with phase drift.

**INPUTS:**

in                = input data in matlab structure format.  
totalDrift       = total amount of phase drift (in degrees) to add over the whole scan.  
noise            = the standard deviation of noise to add to the phase drift function.

#### **4.48. op\_megapressproc.m**

**USAGE:**

```
[out1_diff,out1_sum,outw]=op_megapressproc(filestring,avgAlignDomain,alignSS);
```

**DESCRIPTION:**

Processing script for Siemens MEGA-PRESS MRS data. Includes combination of receiver channels, removal of bad averages, frequency drift correction, manual alignment of edit-on and edit-off spectra, and leftshifting.

**INPUTS:**

filestring = String variable for the name of the directory containing the water suppressed .dat file. Water unsuppressed .dat file should be contained in [filestring '\_w/'];

avgAlignDomain = (Optional) Perform the spectral registration (drift correction) using the full spectrum ('t'), or only a limited frequency range ('f'). Default is 'f'.

alignSS = 0 - Do not align the edit-on and edit-off subspectra. 2 - Perform manual alignment of edit-on and edit-off subspectra.

**OUTPUTS:**

out1\_diff = Fully processed difference spectrum.

out1\_sum = Fully processed sum spectrum.

outw = Fully processed water unsuppressed spectrum.

**4.49. op\_phaseAlignAverages.m****USAGE:**

[out,phs]=op\_phaseAlignAverages(in,Npts,avg,weighting)

**DESCRIPTION:**

Perform time-domain spectral registration using only phase adjustment (no frequency adjustment). This is rarely used.

**INPUTS:**

in = Input data structure.

Npts = Number of points in time domain to use for alignment.

avg = Align averages to the average of the averages ('y'), or the first average in the series ('n');

weighting = (Optional) Apply less weight to the later points of the fid?

**4.50. op\_phaseAlignAverages\_fd.m****USAGE:**

[out,phs]=op\_phaseAlignAverages\_fd(in,minppm,maxppm,Npts,avg,weighting)

**DESCRIPTION:**

Perform time-domain spectral registration using a limited range of frequencies and using only phase adjustment (no frequency adjustment). This is rarely used.

**INPUTS:**

in = Input data structure.

minppm = Minimum of frequency range (ppm).

maxppm = Maximum of frequency range (ppm).

Npts = Number of points in time domain to use for alignment.

avg = Align averages to the average of the averages ('y'), or the first average in the series ('n');

weighting = (Optional) Apply less weight to the later points of the fid?

#### 4.51. op\_pressproc.m

**USAGE:**

```
[out,out_w,out_noproc,out_w_noproc]=op_pressproc(filestring,aaDomain,tmaxin,iterin);
```

**DESCRIPTION:**

Processing script for Siemens PRESS MRS data. Includes combination of receiver channels, removal of bad averages, frequency drift correction, and leftshifting.

**INPUTS:**

**filestring:** String variable for the name of the directory containing the water suppressed .dat file. Water unsuppressed .dat file should be contained in [filestring '\_w/'];

**aaDomain:** (Optional) Perform the spectral registration (drift correction) using the full spectrum ('t'), or only a limited frequency range ('f'). Default is 'f'.

**tmaxin:** (Optional). Duration (in sec.) of the time domain signal used in the spectral registration (drift correction). Default is 0.2 sec.

**iterin:** (Optional). Maximum number of allowed iterations for the spectral registration to converge. Default is 20.

**OUTPUTS:**

**out:** Fully processed, water suppressed output spectrum.

**out\_w:** Fully processed, water unsuppressed output spectrum.

**out\_noproc:** Water suppressed output spectrum without pre-processing (No bad-averages removal, no frequency drift correction).

**out\_w\_noproc:** Water unsuppressed output spectrum without pre-processing.

#### 4.52. op\_readlcmraw.m

**USAGE:**

```
out=readlcmraw(filename,type);
```

**DESCRIPTION:**

Reads LCModel .RAW format into the FID-A data structure format in MATLAB.

**INPUTS:**

**filename** = filename of LCModel raw file.

**type** = type of LCModel raw file:

- 'rda' - .raw file generated from Siemens RDA file
- 'dat' - .raw file generated by FID-A using op\_writelcm.
- 'sim' - .raw file generated from FID-A simulated data.
- 'raw' - not sure about this one.

#### 4.53. op\_readlcmraw\_basis.m

**USAGE:**

```
out=readlcmraw_basis(filename);
```

**DESCRIPTION:**

Reads LCModel .raw model spectrum file into the FID-A data structure format in

MATLAB.

INPUTS:

filename = filename of LCModel raw file.

#### 4.54. op\_rmNworstaverages.m

USAGE:

```
[out,metric,badAverages]=op_rmNworstaverages(in,n);
```

DESCRIPTION:

Removes motion corrupted averages from a dataset containing multiple averages. The N most badly motion corrupted averages are discarded.

INPUTS:

in = input data in matlab structure format  
n = number of bad averages to remove

#### 4.55. op\_rmbadaverages.m

USAGE:

```
[out,metric,badAverages]=op_rmbadaverages(in,nsd,domain);
```

DESCRIPTION:

Removes motion corrupted averages from a dataset containing multiple averages. Bad averages are identified by calculating a 'likeness' metric for each average. This is done by subtracting each average from the average of the averages, and then calculating the root mean squared of this difference spectrum. Averages whose likeness metrics are greater than 'nsd' above the mean are discarded.

INPUTS:

in = input data in matlab structure format  
nsd = number of standard deviations to use a rejection threshold  
domain = domain in which to perform calculations ('t' or 'f')

#### 4.56. op\_rmworstaverage.m

USAGE:

```
[out,metric,badAverages]=op_rmworstaverage(in,n);
```

DESCRIPTION:

Removes motion corrupted averages from a dataset containing multiple averages. The most badly motion corrupted average is discarded.

INPUTS:

in = input data in matlab structure format

#### 4.57. op\_specialproc.m

USAGE:

```
[out,out_w,out_noproc,out_w_noproc]=op_specialproc(filestring,aaDomain,tmaxin,
```

```
iterin);
```

#### DESCRIPTION:

Processing script for Siemens SPECIAL MRS data. Includes combination of receiver channels, removal of bad averages, frequency drift correction, and leftshifting.

#### INPUTS:

filestring: String variable for the name of the directory containing the water suppressed .dat file. Water unsuppressed .dat file should be contained in [filestring '\_w/'];

aaDomain: (Optional) Perform the spectral registration (drift correction) using the full spectrum ('t'), or only a limited frequency range ('f'). Default is 'f'.

tmaxin: (Optional). Duration (in sec.) of the time domain signal used in the spectral registration (drift correction). Default is 0.2 sec.

iterin: (Optional). Maximum number of allowed iterations for the spectral registration to converge. Default is 20.

#### OUTPUTS:

out: Fully processed, water suppressed output spectrum.

out\_w: Fully processed, water unsuppressed output spectrum.

out\_noproc: Water suppressed output spectrum without pre-processing (No bad-averages removal, no frequency drift correction).

out\_w\_noproc: Water unsuppressed output spectrum without pre-processing.

### 4.58. op\_specialproc\_fmrs.m

#### USAGE:

```
[out_stimOFF,out_stimON,out_w]=op_specialproc_fmrs(filestring,blockDesign,leadingAvgstoRmv);
```

#### DESCRIPTION:

Processing script for functional MRS data acquired using the SPECIAL MRS sequence. Includes combination of receiver channels, removal of bad averages, frequency drift correction, and leftshifting. Also includes prior knowledge of the stimulation paradigm (given by the 'blockDesign' vector) and returns the averaged stimulus OFF and stimulus ON spectra.

#### INPUTS:

filestring: String variable for the name of the directory containing the water suppressed .dat file. Water unsuppressed .dat file should be contained in [filestring '\_w/'];

blockDesign: This is a vector of positive and negative even integers that make up the ON/OFF block design. Each integer represents the number of sequential averages in a block. Positive integers refer to ON blocks, and negative integers refer to OFF blocks. For example, for a block design consisting of 30 OFF averages followed by 20 ON averages followed by 10 OFF averages, the blockDesign vector would be: [-30 20 -10];

leadingAvgstoRmv: The number of averages to omit from the beginning of each block. This is done to account for a lag in the neurochemical response to stimulus. Must be an even integer. (optional. Default=0);

#### OUTPUTS:



out\_stimOFF: Fully processed water suppressed spectrum from the sum of the stimulus OFF periods.  
out\_stimON: Fully processed water suppressed spectrum from the sum of the stimulus ON periods.  
out\_w: Fully processed, water unsuppressed output spectrum.

#### **4.59. op\_specialproc\_fmrs\_slidingWindow.m**

**USAGE:**

```
[out1,out_w] = op_specialproc_fmrs_slidingWindow(filestring>windowSize);
```

**DESCRIPTION:**

Processing script for functional MRS data acquired using the SPECIAL MRS sequence. Includes combination of receiver channels, removal of bad averages, frequency drift correction, and leftshifting. This code generates a 'sliding window timecourse' of MR spectra by combining the averages within a small window given by the windowSize argument, and then sliding the window by 1 average and combining again. Each summed window is output as an LCModel text file to be analyzed in LCModel. As a result, this function generates many text output files.

**INPUTS:**

filestring: String variable for the name of the directory containing the water suppressed .dat file. Water unsuppressed .dat file should be contained in [filestring '\_w/'];  
windowSize: This is an integer that specifies the number of averages that are stored within the sliding window. It is recommended to choose a window size that is divisible by the number of phase cycles so that the window does not contain any partial phase cycles.

**OUTPUTS:**

out1: The first sliding window spectrum.  
out\_w: Fully processed, water unsuppressed output spectrum.

#### **4.60. op\_subtractScans.m**

**USAGE:**

```
out=op_subtractScans(in1,in2);
```

**DESCRIPTION:**

Subtract input 2 from input 1;

**INPUTS:**

in1 = 1st input data in matlab structure format.  
in2 = 2nd input data in matlab structure format.

#### **4.61. op\_takeaverages.m**

**USAGE:**

```
out=op_takeaverages(in,index);
```

**DESCRIPTION:**

Extract the averages with indices corresponding to the 'index' input array.

**INPUTS:**

in = input data in matlab structure format.  
index = vector indicating the indices of the averages you would like to extract.

#### **4.62. op\_takesubspec.m**

**USAGE:**

out=op\_takesubspec(in,index);

**DESCRIPTION:**

Extract the subspectra with indices corresponding to the 'index' input array.

**INPUTS:**

in = input data in matlab structure format.  
index = vector indicating the indices of the subspectra you would like to extract.

#### **4.63. op\_unfilter.m**

**USAGE:**

out=op\_unfilter(in,lb);

**DESCRIPTION:**

Multiply the fid by an inverted exponential decay function to undo the effects of filtering.

**INPUTS:**

in = input data in matlab structure format.  
lb = line narrowing factor (Hz).

#### **4.64. op\_writejmrui.m**

**USAGE:**

RF=op\_writejmrui(in,outfile);

**DESCRIPTION:**

Takes MRS data in matlab structure format and writes it to a text file that can be read by jMRUI.

**INPUTS:**

in = input data in matlab structure format.  
outfile = Desired filename of output text file.

#### **4.65. op\_writelcm.m**

**USAGE:**

RF=op\_writelcm(in,outfile,te);

**DESCRIPTION:**

Takes MRS data in matlab structure format and writes it to a text file that

can be read by LCModel.

**INPUTS:**

in               = input data in matlab structure format.  
outfile         = Desired filename of output text file.  
te               = Echo time of acquisition (in ms).

#### **4.66. op\_zeropad.m**

**USAGE:**

out=op\_zeropad(in,zpFactor);

**DESCRIPTION:**

Apply zeropadding (a.k.a. zero-filling) to MRS data.

**INPUTS:**

in               = input data in matlab structure format.  
zpFactor        = the factor by which the number of points in the fid will be increased. ie. if zpFactor =2, then the number of zeros added to the end of the fid will be equal to the number of points in the original spectrum.

#### **4.67. op\_zerotrim.m**

**USAGE:**

out=op\_zerotrim(in,numPointsToTrim);

**DESCRIPTION:**

Remove zeros (or even non-zero data points) from the end of the fid.

**INPUTS:**

in               = input data in matlab structure format.  
numPointsToTrim = The number of points to trim from the end of the fid.

#### **4.68. readlcmcoord.m**

**USAGE:**

out = readlcmcoord(filename,metab)

**DESCRIPTION:**

Reads a LCModel .coord file and extracts the desired part.

**INPUTS:**

filename        = filename of the LCModel .coord file.  
part            = Which metabolite fit to extract from the .coord file - The abbreviated metabolite name should be given (ie. 'Cr', 'PCr', 'Glu', 'GABA', etc.)

#### **4.69. readlcmcoord\_getBackground.m**

**USAGE:**

out = readlcmcoord\_getBackground(filename,part)

DESCRIPTION:

Reads a LCModel .coord file and extracts the desired part.

INPUTS:

filename = filename of the LCModel .coord file.  
part = Which part of the .coord file to extract - 'bg' extracts the  
LCModel baseline signal, 'sp' extracts the spectrum, and 'fit'  
extracts the fit.

#### **4.70. readlcmtab.m**

USAGE:

out = readlcmtab(filename)

DESCRIPTION:

Reads a LCModel .table output file and stores the metabolite concentrations into a matlab structure array.

INPUTS:

filename = filename of the LCModel .table file.